STAY SAFE N

Testing Decedents for COVID-19

04/01/2022

People who have signs and symptoms consistent with COVID-19 may die for a variety of reasons prior to being tested for SARS-CoV-2, the virus that causes COVID-19. The Minnesota Department of Health (MDH) has worked with medical examiners, coroners, and morticians to provide post-mortem SARS-CoV-2 testing supplies by request, as outlined by CDC: Collection and Submission of Postmortem Specimens from Deceased Persons with Confirmed or Suspected COVID-19 (www.cdc.gov/coronavirus/2019-ncov/hcp/guidance-postmortem-specimens.html).

This guidance is intended for morticians, funeral homes, and crematories and provides instructions on how to request supplies from MDH and administer COVID-19 tests to decedents. Questions about this guidance can be directed to Health.Mortsci@state.mn.us.

Requesting COVID-19 specimen collection supplies for a decedent

If a decedent is suspected to have died of COVID-19, a funeral home can request COVID-19 specimen collection supplies from MDH. Decedents tested for COVID-19 should not be embalmed or have their nasopharyngeal passages disinfected prior to specimen collection.

- Funeral homes can request COVID-19 test supplies for decedents by contacting <u>Leslie.kollmann@state.mn.us</u> at MDH. MDH will provide specimen collection supplies and guidance on specimen collection and will coordinate return shipping or specimen pick-up.
- Permission from the family must be obtained prior to requesting specimen collection supplies.
- Have the following information ready prior to requesting specimen collection supplies:
 - Name of the decedent.
 - Date of death.
 - Point of contact for the family, including name, relationship to the decedent, and phone number.
 - Point of contact for the funeral home, including name, and phone number.

If MDH is notified of a possible COVID-19 death that is without a recent positive test, the funeral home may receive a phone call from MDH regarding obtaining a sample from the decedent.

Collecting a COVID-19 respiratory specimen

Morticians who are collecting specimens for COVID-19 testing should use appropriate personal protective equipment (double gloves, fluid resistant gown; plastic face shield or face mask; and goggles) to minimize

TESTING DECEDENTS FOR COVID-19

exposure to infectious bodily fluids. Samples must be collected by a licensed mortician; MDH staff is not able to collect the sample from the decedent.


Specimens may be collected at any time post-mortem, regardless of the extent of body decomposition, as long as they are collected prior to embalming or disinfection of the head or nasopharyngeal passages.

Please collect the specimens for COVID-19 testing before:

- Embalming the body.
- Disinfecting the decedent's head or nasopharyngeal passages.

Collecting a postmortem nasopharyngeal (NP) swab

- 1. Put on appropriate personal protective equipment (gloves, gown, face shield or face mask).
- 2. Label the tube with the decedent's name and date of birth.
- 3. There is one long swab. Choose a nostril and insert parallel to the palate (see image below). You should stop when you feel resistance usually about halfway between the nose and the start of the ear.



- 4. Let the swab sit for a few seconds.
- 5. Remove the swab from the nostril and place it into the tube with the pink or clear liquid (viral transport medium).
- 6. You will notice that the swab has a small perforation about halfway up the stick. Break off the top portion of the swabs at this perforation.
- 7. Close the lid of the tube. It is a tight fit, so you will need to screw the lid on with a bit of pressure.
- 8. Label the tube with the name, date of birth of the person being swabbed and the current date.
- 9. Place the tube into the bag.
- 10. Remove personal protective equipment and perform good hand hygiene.
- 11. Return specimen to MDH as directed.



Minnesota Department of Health | health.mn.gov | 651-201-5000 625 Robert Street North PO Box 64975, St. Paul, MN 55164-0975

Contact <u>health.communications@state.mn.us</u> to request an alternate format.